

ISOLATION OF ACTIN**Number:****Date:****Amount of acetone powder: _____ g**

Check:

1. Extraction with buffer A

	2 l	4 l	
a. Buffer A:			_____
2 mM imidazole-HCl	0.485 g	0.97 g	
0.2 mM Na ₂ ATP	0.22 g	0.44 g	
0.5 mM β-mercaptoethanol	0.078 g	0.156 g	
0.2 mM CaCl ₂	0.059 g	0.118 g	
0.005% NaN ₃	0.1 g	0.2 g	
pH=8.0 (25°C)			
b. Volume of buffer A: 20 ml/g-acetone powder			
c. Extraction: stirring for 30 min, 0-0.5°C			

2. Filtration

- a. Filtration through several layers of sterile cheesecloth and/or
[b. low-speed centrifugation at 500-1,000 g for 10-20 min.]

3. Reextraction of residue

- a. Extraction in 20 ml/g buffer A; 10 min, 0-0.5°C
b. Combine extracts

4. Centrifugation

- a. 10,000-20,000 g, 1 hr, 4°C (Typically, 20,000 rpm, Beckman 55.2 rotor, 30 min)
b. Pipet off (and use) supernatant

5. Polymerization

- a. Concentrations brought to:
- | | |
|------------------|-------|
| KCl (MW: 74.6) | 50 mM |
| Mg ²⁺ | 2mM |
| ATP (MW: 551.1) | 1 mM |
- b. Assemble for 2 hrs, 4°C (on ice)

6. High salt wash (Tropomyosin removal)

- a. Add: solid KCl to 0.6 M
b. Stir gently for 0.5 hr, 4°C (30 minutes from the complete dissolution of KCl)

7. Sedimentation of F-actin

- a. Centrifugation (in 30 ml centrif. tubes), 80,000 g, 3 hrs, 4°C (Typically, 45,000 rpm, Beckman 55.2 rotor, 2 hours)
[b. homogenization of pellet into 150 ml wash buffer (buffer A + 0.6 M KCl, 2 mM MgCl₂, 1 mM ATP)
c. Resediment F-actin 80,000 g, 3 hrs, 4°C]

d. Rinse F-actin pellet thoroughly with buffer A

8. Depolymerization

- a. Allow F-actin pellet to stand in centrifuge tube in ~1 ml/g-acetone powder cold buffer A for 1 hr.
- b. Resuspend F-actin pellet by gentle homogenization in cold buffer A (total A-buffer volume = 3 ml/g-acetone powder)
- c. Dialyze against buffer A, 4°C; change buffer A frequently (e.g., 1st, 3rd, 5th hour, plus overnight)
[for large acetone-powder quantities: divide homogenate into equal 6-ml aliquots and place them into dialysis bags 1/4 inch in diameter and 12 inches long; tie them to magnetic stirring rod in a 1-liter graduated cylinder; dialyze for 6 hrs, while stirring]

9. Clarification of G-actin

- a. Centrifugation of dialyzed actin at 80,000 g, 3 hrs, 4°C
(Typically, 45,000 rpm, Beckman 55.2 rotor, 2 hours)

10. Polymerization

- a. Add:

KCl	50 mM
MgCl ₂	1 mM
ATP	1 mM
NaN ₃	0.02 %
- b. Store at 4°C (on ice; ensure frequent, reliable change of ice)

11. Protein concentration: $[actin](mg/ml) = \frac{OD_{280}}{1.1} \times (dilution)$

Actin cc=

*Note: a) measure against buffer A. b) start with undiluted actin; if OD is in the range of 2, then start diluting with buffer A. c) attach spect printout here.

12. Electrophoresis:

Sample preparation:

Gel loading:

1. 2. 3. 4. 5. 6. 7. 8. 9. 10.

13. Notes/Modifications: