

Preparation of native titin from rabbit *longissimus dorsi*(Based on Soteriou, A., Gamage, M, and Trinick, J. *J. Cell Sci.* **104**, 119-123, 1993)**Number:****I. Solutions**

<i>A. Homogenization buffer</i>	500 ml	1000 ml	2000 ml
50 mM KCl	1.864 g or 10 ml stock (2.5 M)	3.728 g	7.456 g
5 mM EGTA	5 ml stock (0.5 M)	10 ml (0.5 M)	20 ml (0.5 M)
1 mM NaHCO ₃ , pH 7.0	1 ml stock (0.5 M)	2 ml (0.5 M)	4 ml (0.5 M)
0.01% NaN ₃	0.5 ml stock (10%)	1 ml (10%)	0.2 g or 2 ml (10%)

Add fresh:

1 mM DTT	0.5 ml stock (1 M)	1 ml (1M)	2 ml stock (1 M)
1 mM PMSF	2.5 ml stock (0.2 M in EtOH)	5 ml (0.2 M)	10 ml stock
20 µg/ml Trypsin inhibitor	10 mg	20 mg	40 mg
40 µg/ml Leupeptin	1 ml stock (20 mg/ml) =10 aliq.	(10 aliq.)	
20 µM E-64	1 ml stock (10 mM) =10 aliq.	(10 aliq.)	

*Note: prepare 2000 ml's. Take out 500, and complete with E-64 and Leupeptin. To the remaining 1500 ml add 6 aliquots Leupeptin (8 µg/ml) and 6 aliquots E-64 (4 µM).

<i>B. Extraction Buffer</i>	250 ml
10 mM Imidazole-HCl, pH 7.0	2.5 ml stock (1 M stock solution)
0.9 M KCl	16.776 g
2 mM MgCl ₂	0.5 ml stock (1 M)
2 mM EGTA	1 ml stock (0.5 M)
0.01% NaN ₃	0.25 ml stock (10%)

Add fresh:

1 mM DTT	250 µl stock (1 M)
1 mM PMSF	2 ml
40 µg/ml leupeptin	500 µl stock (20 mg/ml) -> 5 aliquots
30 µM E64 (MW=357.4)	500 µl stock (10 mM) -> 5 aliquots
40 µg/ml Trypsin inhibitor	10 mg

<i>C. Precipitation buffer (fresh)</i>	3000 ml	1500 ml
0.1 mM NaHCO ₃	0.6 ml stock (0.5 M)	0.3 ml (0.5 M)
0.3 mM DTT	0.9 ml stock (1 M)	0.45 ml (1 M)
0.1 mM EGTA	0.6 ml stock (0.5 M)	0.3 ml (0.5 M)
2 µg/ml Leupeptin	0.3 ml stock (20 mg/ml) -> 3 aliquots	1.5 aliq.
1 µM E-64	0.3 ml stock (10 mM) -> 3 aliquots	1.5 aliq.

<i>D. Chromatography Buffer</i>	2000 ml	1000 ml	500 ml	400 ml
30 mM K-phosphate, pH 7				
0.6 M KCl	89.472 g	44.736 g	22.368 g	17.894 g

Preparation:

- Prepare 30 mM KH₂PO₄ and K₂HPO₄ buffers separately:

	2000 ml	1000 ml
a. 30 mM KH ₂ PO ₄	8.165 g	4.0825 g
b. 30 mM K ₂ HPO ₄	10.451 g	5.2255 g
- Mix the two so that a pH of 7 is reached. (Lot of dibasic will be needed!).
- Add KCl to the pH-7 phosphate buffer.

Add fresh:

<i>1. For elution:</i>	500 ml
0.1 mM EGTA	0.1 ml (0.5 M)
0.3 mM DTT	0.15 ml stock (1 M)
2 µg/ml Leupeptin	50 µl stock (20 mg/ml) -> 0.5 aliq.
1 µM E-64	50 µl stock(10 mM) -> 0.5 aliq.
0.01% NaN ₃	0.5 ml stock (10%)
0.05% Tween-20	2.5 ml stock (10%)

2. For solubilization of titin pellet:

	10 ml
1 mM EGTA	20 µl stock (0.5 M)
1 mM DTT	10 µl stock (1 M)
20 µg/ml Leupeptin	10 µl stock (20 mg/ml)

0.01% NaN₃ 10 μ l
 0.05% Tween-20 50 μ l

E. Distilled water: 4000 ml Chill overnight!

F. Grinder & Homogenizer soaking solution 2000 ml
 1 mM EDTA (pH 7.4) 4 ml (0.5 M stock)
 0.1 mM PMSF (add fresh) 1 ml (0.2 M stock)

G. Inhibitor Cocktail for titin storage (fresh)
 for 0.1 % NaN₃ 400 μ l stock (10 %)
 for 40 μ g/ml Leupeptin 80 μ l stock (20 mg/ml)
 for 20 μ M E-64 80 μ l stock (10 mM)

(add 14 μ l to each 1-ml fraction)

II. Preparation steps:

Check:

1. Very quickly excise rabbit back muscle.

*Rabbit anaesthesia:

Rabbit weight:	Volume of added anaesthetics
35 mg/kg Ketamin (Calypsol)	(3 ml typ)
7.5 mg/kg Xylozine	(0.5 ml typ)
0.02 mg/kg Atropin	(0 typ)

In deep anaesthesia of the rabbit, incise carotids and bleed rabbit. Soak fur on back of rabbit with cold water. Make incision along spine. Expose back muscles, remove fascia, peel off longissimus dorsi. Put muscle immediately in pre-chilled, pre-washed/soaked electric grinder.

2. Grind and weigh muscle quickly. Weight of muscle (typ. 111 g): _____

Homogenization buffer volume (3 x g-muscle) (typ. 333 ml): _____

Homogenize in 1-l beaker, on ice (ULTRA-TURRAX). _____

3. Wash-centrifuge myofibrils in Homogenization buffer 3 times (2000g) _____

(three tubes in Sorvall GSA rotor, 4500 rpm, 10 min, 0 °C)

Buffer volume (3 x g-muscle) [typ. 111 ml/each of three tubes]: _____

4. Suspend final pellet in Extraction solution (2 x g-muscle) _____

(Do not remove pellet from centrifuge tube.)

Buffer volume [typ. 74 ml/each of three tubes]: _____

5. Extract for 5 minutes while stirring. Do procedure in Centrifuge tube! _____

6. Centrifuge 20,000g, 30 min. (Sorvall GSA rotor, 12,000 rpm, 30 min, 0 °C) _____

7. Supernatant diluted 3 times (to achieve 0.2 M ionic strength) _____

Supernatant volume (typ ~90 ml): _____

Volume of added Precipitation solution (typ ~180 ml): _____

Total volume (typ ~270 ml): _____

*Note: optionally, filter supernatant through gauze to get rid of fat particles.

8. Incubate on ice for 1 hour _____

9. Centrifuge, 20,000g, 30 min (sediment myosin) _____

(Sorvall GSA rotor, 12,000 rpm, 30 min, 0 °C)

10. Supernatant diluted 5 times (final ionic strength 0.05 M) _____

Supernatant volume (typ ~245 ml): _____

Volume of added Precipitation solution (typ ~980 ml): _____

Total volume (typ ~1225 ml): _____

11. Incubate on ice for 40 minutes _____

12. Centrifuge, 10,000g, 30 min (Sorvall GSA rotor, 8000 rpm, 30 min, 0 °C) _____

13. Titin pellet solubilized in Chromatography buffer

Volume of Chromatography buffer (typ ~700 μ l/tube+2x700 μ l wash):

Use plastic transfer pipet to softly dislodge pellet. Pellet is on side of tube(!)

Add as much chromatography buffer so that OD is about 15.

14. Centrifuge, 25,000g, 30 min (Sorvall UC, 17,000 rpm, 30 min, 0°C) _____

Supernatant volume (typ ~4 ml):

15. Supernatant loaded onto Sepharose CL-2B column. _____

Column: Bio-Rad 1x120 cm, max. vol. 103 ml

Chromatography Conditions:

a. Sample loaded onto column (typ in two 2.0-ml turns, labelled as #a and #b):

-Volume (typ 2.0): #a: #b:

-OD (typ ~18/ml):

b. Flow rate(Typ 0.2 ml/min; 80 setting on Gilson pump using narrowest tube):

c. Fraction size (Typ 1ml):

d. Elution/fraction collection:

-Void volume (Typ 20 ml): #a: #b:

-Paper speed: (Typ 0.2 mm/min):

-AUFS (Typ 0.2):

e. Elution profile:

Abs₂₈₀:

-UV-Monitor: attach records

-Spectrophotometer (260 nm, 280 nm): attach records

16. SDS-PAGE:

Gel type: Fairbanks (Typ 2.5-12%). Note: FBX denotes Fairbanks solubilization sol.

Sample preparation:

a. Standard: Rabbit LD

10 μ l LD

+90 μ l 1xFBX (incubate on 65 °C for 45 s)

Load 3-5 μ l on gel.

b. Titin fractions:

-10 μ l sample + 5 μ l 3xFBX + 35 μ l 1xFBX (incubate on 65 °C for 45 s).

-10 μ l of the above mixture + 90 μ l 1xFBX (incubate on 65 °C for 45 s).

Load 20 μ l on gel.

[**17.** Optionally, pool titin peak. Use in *in vitro* motility assay]

18. Store titin:

To the chromatography fractions, add:

-0.1% NaN₃

-40 μ g/ml Leupeptin

-20 μ M E-64

Then, a. Quick-freeze in liquid nitrogen, store at -80 °C

or, b. Store some on ice for a few days.

III. Notes, Remarks, Modifications: