

Preparation of myofibrils from glycerinated psoas muscle

Knight, P.J. & Trinick, J.A., (1982) Methods Enzymol. 85: 9

Number:

Buffers:

Rigor buffer:

0.1 M KCl
2 mM MgCl₂
1 mM EGTA
0.5 mM DTT (fresh)
10 mM K-phosphate
pH 7 (at 0°C)

Vol: 1000 ml

7.456 g
2 ml stock (1 M)
2 ml stok (0.5 M)
0.5 ml stock (1 M)
prepare from concentrated stock

Na-Rigor buffer:

0.1 M NaCl
2 mM MgCl₂
1 mM EGTA
0.5 mM DTT (fresh)
10 mM K-phosphate
pH 7 (at 0°C)

Vol: 1000 ml

5.844 g
2 ml stock (1 M)
2 ml stok (0.5 M)
0.5 ml stock (1 M)
prepare from concentrated stock

Steps:

1. Tease muscle into fine fiber bundles
2. Soak bundles in rigor buffer for 1 hr (Rotorack in 15 ml Falcon tube)
Approximate muscle volume= ml
Volume of added rigor buffer
3. Homogenize in 10 vol. rigor buffer (Ultra-Turrax on ice,
in 50 ml Falcon tube)
4. Sediment at 2000g for 5 min. (3000 rpm, Janetzky large rotor, 4 °C)
in 15 ml Falcon tube
5. Wash pellet in 10 vol. rigor buffer
6. Repeat steps 4-5 two to three times

Measuring protein concentration:

1. Centrifuge an aliquot of myofibrils and resuspend the pellet in Na-rigor buffer.
2. Dissolve suspension in warm 1% SDS.
3. Measure absorbance at 280 nm. Extinction coefficient: 0.7 ml/mg-cm.

Storage:

1. Sediment myofibrils from suspension by centrifugation.
2. Resuspend pellet with rigor buffer containing 50% glycerol.
3. Store at -20°C.