

WORKSHEET**Papain Mg-S1 preparation****Number:****Date:****Number of myosin preparation:****Solutions, buffers:****Check:**

1. BED	Volume: 100 ml	_____
0.1 mM NaHCO ₃	0.02 ml stock (0.5M)	
0.1 mM EGTA	0.04 ml stock (0.25M) [20 μ l 0.5 M]	
1 mM DTT	100 μ l (fresh 1 M stock)	
2. PMB	Volume 25 ml	_____
25 mM imidazole-HCl (pH 7.4)	0.625 ml stock (1 M)	
100 mM NaCl	1.25 ml stock (2 M) [0.1461 g]	
5 mM MgCl ₂	62.5 μ l stock (2 M) [125 μ l 1M]	
1 mM DTT	25 μ l stock (1 M)	
3. E64 stock	Volume: 1 ml	_____
E64	1 mg	
DMSO	to 1 ml	
*Note: alternatively, use water-based stock. Record stock concentration.		
4. Papain	Volume: 2 ml	_____
5 mg/ml papain	0.01 g	
Papain Base Buffer	to 2 ml	
5. Inhibitor buffer	Volume: 10 ml	_____
5 mM MgCl ₂	25 μ l stock (2 M) [50 μ l 1M]	
5 μ g/ml E64	50 μ l (stock)	
BED	to 10 ml	
6. 1.2 M KCl + 20 mM DTT/BED	Volume: 10 ml	
1.2 M KCl	0.9 g	
20 mM DTT	200 μ l stock (1 M)	
BED	to 10 ml	
7. Papain Base Buffer (pH 7.4 with NaOH)	Volume: 100 ml	
5 mM cysteine	0.06 g	
2 mM EDTA	0.058 g [0.4 ml stock (0.5M)]	

Steps:**Check**

1. Add 9 volumes of BED to 20 mg stock myosin in a centrifuge tube (Falcon tube, with volume marks), and mix to precipitate the filaments. Incubate on ice for >10 min.

2. Centrifuge at low speed in a swinging bucket rotor. _____
(Typically, 1,600 rpm, old IEC, 1 hour)
3. Dissolve the pellet in equal volume of 1.2 M KCL+20 mM DTT _____
in BED.
Pellet volume (estimate from marks on Falcon tube):
4. Add 19x volume of BED, incubate on ice for 10 min, centrifuge _____
(Typically, 1,600 rpm, old IEC, 1 hour)
5. Resuspend pellet in PMB to a final concentration of 10-12 mg/ml, _____
incubate for 10 min at 25°C.
Pellet volume:
Volume of added buffer:
Myosin concentration ($A_{280} = 0.53 \text{ cm}^2/\text{mg}$):
6. Add papain to 12.5 $\mu\text{g}/\text{ml}$, mix gently, incubate 7.5-10 min at 25°C. _____
Volume of papain stock added:
Digestion time:
7. Add greater than equal volume of inhibitor solution to stop reaction _____
Volume of inhibitor added:
8. Incubate on ice for 1 hr. _____
9. High-speed centrifugation. _____
(Typically, 40,000 rpm, Beckman 50Ti, 30 min, 0 °C)
10. Store supernatant (typically 1.1 mg/ml) on ice. _____
Use in the motility assay for 3-5 days. $A_{280}=0.81 \text{ cm}^2/\text{mg}$.

Concentration: $[S - 1](\text{mg} / \text{ml}) = \frac{OD_{280}}{0.81} \times (\text{dilution})$

*Note: attach spect printout here.

SDS-PAGE:

Type of gel:
Percentage:
Sample preparation:

Standard:
Loading: lanes

1. 2. 3. 4. 5. 6. 7. 8. 9. 10.

Notes/Modifications: