

WORKSHEET HMM preparation

Number:

Myosin preparation used:

Solutions, buffers:

Check:

| | | |
|--------------------------------------|-------------------------|-------|
| 1. BED | Volume: 100 ml | _____ |
| 0.1 mM NaHCO ₃ | 0.02 ml stock (0.5M) | |
| 0.1 mM EGTA | 0.04 ml stock (0.25M) | |
| 1 mM DTT (add fresh from 1 M stock) | 100 µl | |
| 2. 2xCHB | Volume 25 ml | _____ |
| 20 mM imidazole-HCl (pH 7.0) | 0.5 ml stock (1M) | |
| 1 M KCl | 10 ml stock (2.5 M) | |
| 4 mM MgCl ₂ | 0.05 ml stock (2 M) | |
| 10 mM DTT (add fresh from 1 M stock) | 0.25 ml | |
| 3. PMSF stock | Volume: 10 ml | _____ |
| 0.2 M PMSF | 0.348 g | |
| Ethanol | to 10 ml | |
| 4. Chymotrypsin | Volume: 2 ml | _____ |
| 5 mg/ml TLCK-chymotrypsin | 0.01 g | |
| 2xCHB | to 2 ml | |
| <i>OR:</i> | Volume: 1 ml | |
| 0.5 mg/ml TLCK-chymotrypsin | 0.5 mg | |
| 2xCHB | to 1 ml | |
| 5. Inhibitor buffer | Volume: 25 ml | _____ |
| 3 mM MgCl ₂ | 10 ml | |
| | 0.075 ml stock (1 M) | |
| | 30 µl (1 M) | |
| 0.1 mM PMSF | 0.0125 ml stock (0.2 M) | |
| BED | 5 µl (0.2 M) | |
| | to 25 ml | |
| | to 10 ml | |

Steps:

- Add 9 volumes of BED to stock myosin in a centrifuge tube, and mix to precipitate the filaments. (Slowly pipet buffer on myosin while vortexing.) Incubate on ice for >10 min.
Volume of myosin: _____
Total volume: _____
- Centrifuge at low speed in a swinging bucket rotor.
(2-4000 rpm, 30 min, 0 °C) _____
- Dissolve the pellet in 2xCHB, and add BED as needed to achieve a final concentration of ~15 mg/ml myosin. (Typically, one-to-one volume ratio)
Volume of pellet (est.): _____
Volume of added 2xCHB _____
- Incubate myosin solution at 25°C for 10 min.
(On water bath. Monitor temperature.) _____
- Add TLCK-treated chymotrypsin to a final concentration of 12.5 µg/ml, gently mix the reaction and incubate 7.5 to 10 min at 25°C.
Volume of added chymotrypsin stock: _____
Digestion time: _____
- Add 9 volumes of ice cold inhibitor buffer to the reaction and mix. _____

(Slowly pipet buffer on protein while vortexing.)
Incubate on ice for 1 hour.

7. Centrifuge the suspension at high speed. _____
[Beckman 50Ti rotor, 40,000 rpm (~100,000 g), 30 min, 0 °C]

8. Store supernatant (typically 0.4-0.7 mg/ml, and >90% pure) on ice. _____
Use in the motility assay for 3-5 days. $A_{280}=0.60 \text{ cm}^2/\text{mg}$.

SDS-PAGE:

Type of gel: Percentage: Sample preparation: Standard:

Loading: 1. 2. 3. 4. 5. 6. 7. 8. 9. 10.

Actin-affinity purification of HMM:

Solutions, buffers:

1. HMM (see above), concentration determined _____
2. F-actin (concentration determined). **Number** of Actin prep: _____
3. ATP stock (300 mM) _____
4. AB buffer (see *In vitro* assay worksheet) _____

Steps:

1. Dilute HMM in AB buffer to a final concentration of about 0.4 mg/ml _____
Volume of HMM:
Volume of added AB:

2. Add F-actin to a final concentration of 1.6 mg/ml. _____
Mix gently but thoroughly.
Volume of F-actin solution:

3. Add ATP stock to a final concentration of 1 mM _____
(or to a final concentration not more than that used in the *in vitro* assay)
Mix gently but thoroughly.
Volume of ATP stock:

4. Incubate mixture on ice for 10 minutes _____

5. Centrifuge _____
[Beckman 50Ti rotor, 45,000 rpm (150,000 g), 30 min, 0 °C]

6. Use supernatant in the *in vitro* assay, diluted typically to 40-60 $\mu\text{g}/\text{ml}$.

Example:

| Component | Volume | Final concentration |
|------------------------------|-----------------|---------------------|
| HMM (stock cc.: 0.67 mg/ml) | 1.2 ml | 0.4 mg/ml |
| F-actin (stock cc.: 5 mg/ml) | 0.64 ml | 1.6 mg/ml |
| ATP (stock cc.: 300 mM) | 6 μl | 1 mM |
| AB buffer | 0.154 ml | |
| Total volume: | 2 ml | |