

Preparation of and *in vitro* motility assay on human cardiac myosin

Specimen number:

Experiment number:

Buffers, solutions:

1. Myosin extraction buffer		volume: 500 ml		
	0.3 M KCl		11.184 g	
	0.09 M KH ₂ PO ₄		6.124 g	
	0.06 M K ₂ HPO ₄		5.225 g (pH 6.8)	
	1 mM MgCl ₂		500 µl stock (1 M)	
Fresh for 5 ml:	0.2 mM ATP		3.3 µl (300 mM stock)	
	1 mM DTT		5 µl (1 M stock)	
	0.23 mM PMSF		5.75 µl (0.2 M stock in EtOH)	
2. AB buffer		Volumes:	100 ml	250 ml
	25 mM imidazole-HCl (pH 7.4)		0.17 g	0.425 g
	25 mM KCl		0.186 g	0.466 g
	4 mM MgCl ₂ (1 M stock)		0.4 ml	1 ml
	1 mM EGTA (0.5 M stock)		0.2 ml	0.5 ml
	1 mM DTT (fresh 1 M stock)		100 µl	250 µl
				0.154 g
3. 0.6 M KCl/AB		Volumes:	50 ml	100 ml
	0.6 M KCl		2.144g	4.287g
	— add AB buffer —			
4. AB/BSA		Volumes:	50 ml	100 ml
	0.5 mg/ml BSA		0.025 g or	0.05 g or
			125 µl stock (200 mg/ml)	250 µl stock
	— add AB buffer —			
5. GOC		Volumes:	2 ml	5 ml
	glucose stock*		40 µl	100 µl
	catalase**		4 µl	10 µl
	glucose oxidase***		16 µl	40 µl
	AB/BSA		1.925 ml	4.815 ml
	β-mercaptoethanol		15 µl	37.5 µl

*glucose stock is 30% (1.665M) in water + NaN₃ stored at 4 °C. Final concentration of glucose is 3 mg/ml.

**catalase is prepared fresh, stock concentration 20 mg/ml in AB buffer.

***glucose oxidase is prepared fresh, stock concentration 25 mg/ml in AB buffer.

6. ATP		Volumes:	1 ml	3 ml	Actual prep:
	1 mM ATP (from 300 mM stock)		3.3 µl	10 µl	
	GOC		0.9967 ml	2.99 ml	

*ATP stock is 300 mM aqueous solution of ATP, pH 6.8, in 50 µl aliquots, quick-frozen in liquid nitrogen, and stored at -20 °C or, preferably, at -80 °C. Final concentration is adjusted as desired.

7. Actin		Volume:	1 ml	5 ml	Actual prep:
	TRITC-Ph-F-actin stock		2 µl	10 µl ml	Prep #:
	ATP		0.998 ml	4.99 ml	

*Note: dilute actin so that a approximate final dilution of 50,000 of the stored F-actin is obtained. Adjust to a concentration that produces an optimal distribution under the fluorescence microscope.

8. Myosin. Preparation depends on final myosin concentration. For the *in vitro* motility assay, 0.1-0.2 mg/ml myosin concentration should be used.

Experiment Steps:

Check:

A. Preparation and storage of muscle sample:

1. Weigh human cardiac muscle sample. Weight = mg. _____
2. Take ~50 mg muscle and place in a 1.6-ml Eppendorf centrifuge tube. Pierce the tube cap to prevent explosion during and after quick freezing. _____
3. In the Eppendorf tube, mince muscle into fine pieces with the tips of a pair of fine scissors. _____
4. Label the tube properly, then quick-freeze the tube in liquid nitrogen. _____
5. Store in -80 °C freezer until preparation of myosin. _____

B. Preparation of myosin:

1. Add 10 ml/g-muscle myosin extraction solution to the muscle (~0.5 ml). Buffer volume:ml. _____
2. Place muscle suspension into glass homogenizer, and homogenize for 30 min. on ice. _____
3. Incubate on Rotorack for 30 min. in cold room. _____
4. Centrifuge at 140.000g, 30 min, 4 °C. _____
5. Dilute supernatant 10-fold with 1 mM DTT in ddH₂O.
Preparation of 1 mM DTT: add 20 µl DTT stock (1 M) to 20 ml water)
Volume of supernatant:ml.
Volume of 1 mM DTT added:ml. _____
6. Incubate on ice for 60 minutes (to precipitate myosin). _____
7. Centrifuge at 20.000g, 20 min, 4 °C. _____
8. Dissolve pellet in ~equal volume of 0.6 M KCl/AB solution.
Volume of 0.6 M KCl/AB added (usually 100 µl):µl. _____
9. Measure concentration of myosin against 0.6 M KCl/AB as background.
$$[\text{myosin}](\text{mg} / \text{ml}) = \frac{OD_{280}}{0.53} \times (\text{dilution})$$

Myosin concentration =mg/ml. _____
10. Dissolve an aliquot of myosin in Fairbanks sample buffer
Fresh Fairbanks sample buffer: 20 mg DTT in 1 ml 3x Fairbanks buffer
2 volumes myosin + 1 volume Fairbanks buffer, 1 min. at 90 °C. Quickfreeze 20 µl aliquots, store in -80 °C freezer. _____
11. Immediately carry out the *in vitro* motility assay. _____

C. Advance preparations for the *in vitro* motility assay:

1. Prepare flow cells (dip-coat 18 mm² coverslips with nitrocellulose (1% in amyl-acetate), and allow to dry overnight. _____
2. Place two Parafilm spacers (4x25 mm) on a precleaned microscope slide, position the 18 mm² coverslip on the top, and melt the Parafilm on a hotplate. _____
3. Store flow cells in a dust-free area. Use within 2-3 days. _____

D. The *in vitro* motility assay:

1. Turn on arc lamp (HBO200) ~10 min prior to experiments. _____
2. Turn on thermostat adjusted to an effective temperature of 37 °C (46 °C on thermostat) _____
3. Place the flow cell obliquely in a 50-ml beaker, and start adding the solutions: _____
 - a. 10 µl myosin, wait 1 min. _____
 - b. 100 µl AB/BSA. _____
 - c. 100 µl actin in ATP/GOC solution _____
4. Oil the objective, and place the flow cell on the microscope stage.
Wait for ~30 s for temperature to equilibrate. _____
5. Turn on the camera and bring actin filaments in focus. Adjust contrast and frame averaging with Argus-20 (usually 4 frames). Start video recording. _____

Temperature:

Tape Number:

Telescope lens setting (default 2x): 1.25x 1.6x 2.0x

Video Recording:

Experiment	Step	Tape Time		
1.	START	h	m	s
		h	m	s
	STOP	h	m	s
2.	START	h	m	s
		h	m	s
	STOP	h	m	s
3.	START	h	m	s
		h	m	s
	STOP	h	m	s
4.	START	h	m	s
		h	m	s
	STOP	h	m	s
5.	START	h	m	s
		h	m	s
	STOP	h	m	s

Measurement of filament velocity

1. Using NIH Image 1.61, grab a sequence of images (10-20 frames).
2. Carry out frame averaging if necessary.
3. Import the "In Vitro Motility Assay 2.3*" macro.
4. Execute the "Preprocess stack" command.
5. Make a rectangular ROI for a selected filament and execute the "AutoMeanFast" command.
6. Repeat step 5 for preferably each filament in the field of view.
7. Export the Results.
8. Import the Results text file into KaleidaGraph, and construct velocity histogram.

Spatial calibration for the 2x telescope lens setting (as of 2000/06/07):

6.374 pixels/ μ m

[For individual-frame measurements:

horizontal scale: 6.591 pixels/micron; vertical scale: 6.157 pixels/micron, pixel aspect ratio: 1.07]

Notes/Remarks/Modifications:

References:

1. Nguyen, T.-T.T., *et al. Circ. Res.* **79**, 222-226, 1996. (*in vitro* motility using human cardiac biopsy samples)
2. Kron, S.J., *et al. Methods Enzymol.* **196**, 399-416, 1991. (*in vitro* motility assay)
3. Margossian, S.S., & Lowey, S. *Methods Enzymol.* **85**, 55, 1982. (skeletal myosin preparation)