

Preparation of Protein-A – anti-DIG-coated latex beads**Type of bead/diameter:****Antibody used: Anti-Digoxigenin, polyclonal, Roche Cat. No. 1 333 089****I. Solutions:***1. Bead suspension:*

3.18 µm diameter carboxylated latex beads (Spherotech), 5 %.

2. Protein-A 5 mg/mlDissolve Protein-A (Sigma) in H₂O at a final concentration of 5 mg/ml.

Prepare 50-µl aliquots.

Quick-freeze aliquots and store at -80 °C.

3. Na-phosphate 100 mlNaH₂PO₄ (monobasic) 0.1 M, pH 5 1.38 g

This monobasic solution will have a pH of ~5.

4. Carbodiimide (EDC or EDAC)

10 mg/ml in Na-phosphate

Prepare at least 1 ml.

5. Borate 100 ml 200 ml

0.1 M Na-borate 3.814 g 7.628 g

Adjust to pH 8.2

6 Triethanolamine 100 ml 200 ml

0.2 M Triethanolamine 3.714 g 7.428 g

Adjust to pH 8.2

7. DMP (dimethylpimelimidate) 10 ml 20 ml

50 mM DMP 0.13 g 0.26 g

Always make fresh.

Dissolve in Triethanolamine

Readjust pH to 8.2.

8. Ethanolamine 10 ml 100 ml 200 ml

50 mM Ethanolamine 31 µl 310 µl 620 µl

Concentration should equal that of DMP.

Dissolve in H₂O

Adjust pH to 8.2.

9. AB buffer Volumes: 100 ml 250 ml 1000 ml

25 mM imidazole-HCl (pH 7.4) 0.17 g 0.425 g 1.7 g

25 mM KCl 0.186 g 0.466 g 1.864 g

4 mM MgCl₂ (1 M stock) 0.4 ml 1 ml 4 ml

1 mM EGTA (0.5 M stock) 0.2 ml 0.5 ml 2 ml

1 mM DTT (add fresh from 1 M stock) 100 µl 250 µl 0.154 g

10. 0.2 M KCl/AB Volumes: 50 ml 100 ml

0.2 M KCl 0.7456 g 1.4914 g

-add AB buffer-

11. Bead blocking solution: Volume: 1 ml

1% BSA 50 µl stock (20%)

0.2 % Tween-20 20 µl stock (10%)

12. PBS 1 liter10 mM K₂HPO₄/KH₂PO₄, pH 7.4

140 mM NaCl 8.18 g

Kellermayer. latex bead coating.

Date:
Number:

Preparation suggestions:

a) Prepare 500 ml 20 mM K-phosphate buffer from mono- and dibasic (pH 7.4). b) Add 8.18 g NaCl to 500 ml 20 mM K-phosphate solution, then bring volume to 1 liter with dH₂O.

II. Steps for protein-A – anti-digoxigenin beads:

A. Dissolution of Anti-Digoxigenin

1. Add 0.4 ml PBS to antibody liophilizate
2. Final antibody concentration is 0.5 mg/ml

B. Preparation of Protein-A-conjugated latex beads (reduce quantities as necessary)

1. Mix: -100 µl Protein-A
 -1 ml Bead suspension
 -1 ml EDC
2. Vortex
3. Rotorack for 2 hours at room temperature.
4. Centrifuge in tabletop centrifuge, 3000 rpm, 4 °C, 4 min.
5. Resuspend pellet in 1 ml 0.2 M KCl AB
6. Repeat steps 4-5 twice. Final bead concentration is 5%.

C. Conjugation of Anti-Digoxigenin antibody to Protein-A beads. (reduce quantities as necessary)

1. Centrifuge 400 µl Protein-A beads in tabletop centrifuge, 3000 rpm, 4 °C, 5 min.
2. Resuspend pellet in 400 µl Borate.
3. Repeat steps 1-2 twice.
4. Add 300 µl ANTI-DIG to the Protein-A bead susedsion.
5. Rotorack for 30 minutes at room temperature.
6. Centrifuge at 3000 rpm, 4 °C, 5 min.
7. Resuspend pellet in 800 µl borate.
8. Centrifuge at 3000 rpm, 4 °C, 5 min.
9. Resuspend pellet in 800 µl Triethanolamine.

D. Croslinking of anti-digoxigenin to Protein-A. (reduce quantities as necessary)

1. Centrifuge ANTI-DIG-conjugated beads at 3000 rpm, 4 °C, 5 min.
2. Resuspend pellet in 800 pl Triethanolamine.
3. Repeat steps 1-2.
4. Centrifuge at 3000 rpm, 4 °C, 5 min.
5. Resuspend pellet in 3 ml DMP. Prevortex pellet, then add DMP while vortexing.
6. Rotorack for 45 minutes at room temperature.
7. Centrifuge at 3000 rpm, 4 °C, 5 min.
8. Resuspend pellet in 800 µl Ethanolamine.
9. Rotorack for 5 minutes at room temperature.
10. Centrifuge at 3000 rpm, 4 °C, 5 min.
11. Resuspend pellet in 800 µl Borate.
12. Centrifuge at 3000 rpm. 4 °C, 5 min.
13. Resuspend pellet in 800 µl 0.2 M KCl-AB.
14. Centriguge at 3000 rpm, 4 °C, 5 min.
15. Repeat steps 13-14.
16. Resuspend pellet in 400 µl 0.2 M KCl-AB. Final bead concentration is ~5 %.

III. Steps for direct bead labeling with anti-digoxigenin:

1. Mix: -100 µl anti-dig
 -0.1 ml Bead suspension
 -0.1 ml EDC
2. Vortex
3. Rotorack for 2 hours at room temperature.
4. Centrifuge in tabletop centrifuge, 3000 rpm, 4 °C, 4 min.
5. Resuspend pellet in 1 ml 0.2 M KCl AB
6. Repeat steps 4-5 twice. Final bead concentration is 5%.