

Preparation of Protein-A-antibody-coated latex beads**Type of bead/diameter:****Antibody used:****I. Solutions:***1. Bead suspension:*

3.18 µm diameter carboxylated latex beads (Spherotech), 5 %.

2. Protein-A 5 mg/mlDissolve Protein-A (Sigma) in H₂O at a final concentration of 5 mg/ml.

Prepare 50-µl aliquots.

Quick-freeze aliquots and store at -80 °C.

3. Na-phosphate 100 mlNaH₂PO₄ (monobasic) 0.1 M, pH 5 1.38 g

This monobasic solution will have a pH of ~5.

4. Carbodiimide (EDC or EDAC)

10 mg/ml in Na-phosphate

Prepare at least 1 ml.

5. Borate 100 ml 200 ml

0.1 M Na-borate 3.814 g 7.628 g

Adjust to pH 8.2

6 Triethanolamine 100 ml 200 ml

0.2 M Triethanolamine 3.714 g 7.428 g

Adjust to pH 8.2

7. DMP (dimethylpimelimidate) 10 ml 20 ml

50 mM DMP 0.13 g 0.26 g

Always make fresh.

Dissolve in Triethanolamine

Readjust pH to 8.2.

8. Ethanolamine 10 ml 100 ml 200 ml

50 mM Ethanolamine 31 µl 310 µl 620 µl

Concentration should equal that of DMP.

Dissolve in H₂O

Adjust pH to 8.2.

9. AB buffer Volumes: 100 ml 250 ml 1000 ml

25 mM imidazole-HCl (pH 7.4) 0.17 g 0.425 g 1.7 g

25 mM KCl 0.186 g 0.466 g 1.864 g

4 mM MgCl₂ (1 M stock) 0.4 ml 1 ml 4 ml

1 mM EGTA (0.5 M stock) 0.2 ml 0.5 ml 2 ml

1 mM DTT (add fresh from 1 M stock) 100 µl 250 µl 0.154 g

10. 0.2 M KCl/AB Volumes: 50 ml 100 ml

0.2 M KCl 0.7456 g 1.4914 g

-add AB buffer-

11. Bead blocking solution: Volume: 1 ml

1% BSA 50 µl stock (20%)

0.2 % Tween-20 20 µl stock (10%)

II. Steps:

A. *Preparation of Protein-A-conjugated latex beads:*

1. Mix: -100 μ l Protein-A
-1 ml Bead suspension
-1 ml EDC
2. Vortex
3. Rotorack for 2 hours at room temperature.
4. Centrifuge in new IEC fuge, 3000 rpm, 4 °C, 4 min.
5. Resuspend pellet in 1 ml 0.2 M KCl AB
6. Repeat steps 4-5 twice. Final bead concentration is 5%.

B. *Conjugation of T12 antibody to Protein-A beads.*

1. Centrifuge 400 μ l Protein-A beads in new IEC, 3000 rpm, 4 °C, 5 min.
2. Resuspend pellet in 400 μ l Borate.
3. Repeat steps 1-2 twice.
4. Add 300 μ l T12 (from Fürst) to the Protein-A bead suspension.
5. Rotorack for 30 minutes at room temperature.
6. Centrifuge at 3000 rpm, 4 °C, 5 min.
7. Resuspend pellet in 800 μ l borate.
8. Centrifuge at 3000 rpm, 4 °C, 5 min.
9. Resuspend pellet in 800 μ l Triethanolamine.

C. *Crosslinking of T12 to Protein-A.*

1. Centrifuge T12-conjugated beads at 3000 rpm, 4 °C, 5 min.
2. Resuspend pellet in 800 μ l Triethanolamine.
3. Repeat steps 1-2.
4. Centrifuge at 3000 rpm, 4 °C, 5 min.
5. Resuspend pellet in 3 ml DMP. Prevortex pellet, then add DMP while vortexing.
6. Rotorack for 45 minutes at room temperature.
7. Centrifuge at 3000 rpm, 4 °C, 5 min.
8. Resuspend pellet in 800 μ l Ethanolamine.
9. Rotorack for 5 minutes at room temperature.
10. Centrifuge at 3000 rpm, 4 °C, 5 min.
11. Resuspend pellet in 800 μ l Borate.
12. Centrifuge at 3000 rpm, 4 °C, 5 min.
13. Resuspend pellet in 800 μ l 0.2 M KCl-AB.
14. Centrifuge at 3000 rpm, 4 °C, 5 min.
15. Repeat steps 13-14.
16. Resuspend pellet in 400 μ l 0.2 M KCl-AB. Final bead concentration is ~5 %.

D. *Pre-blocking T12 beads:*

1. Mix 100 μ l T12-bead suspension and 400 μ l Bead-blocking solution.
2. Rotorack at room temp. for 5 min, vortexing 5s every minute.
3. Centrifuge in new IEC at 3000 rpm, 5 min.
4. Resuspend in 500 μ l 0.2M KCl AB
5. Repeat steps 3-4.

E. *Coating T12 beads with titin*

1. Mix 50 μ l pre-blocked T12 bead suspension, 5-15 μ l titin, and 0.2M KCl AB (to 1 ml).
2. Rotorack at room temp for 30 min.
3. Centrifuge in new IEC at 3000 rpm, 5 min.
4. Resuspend in 1 ml 0.2 M KCl AB. Gently!