

WORKSHEET

Actomyosin in vitro motility assay

(Ref.: Kron, S.J., et al. *Methods Enzymol.* **196**, 399-416, 1991.)

Experiment description:

Experiment Number:

I. Solutions:

1. AB buffer	Volumes:	100 ml	250 ml	1000 ml
25 mM imidazole-HCl (pH 7.4)		0.17 g	0.425 g	1.7 g
25 mM KCl		0.186 g	0.466 g	1.864 g
4 mM MgCl ₂ (1 M stock)		0.4 ml	1 ml	4 ml
1 mM EGTA (0.5 M stock)		0.2 ml	0.5 ml	2 ml
1 mM DTT (add fresh from 1 M stock)		100 μ l	250 μ l	0.154 g
2. 0.2 M KCl/AB	Volumes:	50 ml	100 ml	
0.2 M KCl		0.7456g	1.4914g	
—add AB buffer—				
3. 0.6 M KCl/AB	Volumes:	50 ml	100 ml	
0.6 M KCl		2.237g	4.474g	
—add AB buffer—				
4. AB/BSA	Volumes:	50 ml	100 ml	
0.5 mg/ml BSA		0.025 g	0.05 g	
—add AB buffer—				
5. 0.2 M KCl AB/BSA	Volumes:	50 ml	100 ml	
0.5 mg/ml BSA		0.025 g	0.05 g	
—add 0.2 M KCl AB buffer—				
6. Blocking solution	Volumes:	2 ml	5 ml	
5% BSA		0.5 ml stock*	1.25 ml stock*	
1% Gelatin			100 μ l stock**	250 μ l stock**
0.2% (v/v) Tween-20			40 μ l stock***	100 μ l stock***
—add 0.2 M KCl AB buffer—			1.36 ml	3.4 ml

*BSA stock is 20% (w/v) in 0.2 M KCl AB buffer with 0.1% NaN₃.

**Gelatin stock is 20% (w/v) in 0.2 M KCl AB buffer with 0.1% NaN₃.

***Tween-20 stock is 10% (v/v) in water.

7. GOC	Volumes:	2 ml	5 ml	10 ml
glucose stock*		40 μ l	100 μ l	200 μ l
catalase**		4 μ l	10 μ l	20 μ l
glucose oxidase***		16 μ l	40 μ l	80 μ l
AB/BSA		1.925 ml	4.815 ml	9.63 ml
β -mercaptoethanol		15 μ l	37.5 μ l	75 μ l

*glucose stock is 30% (1.665M) in water stored on ice. Final concentration of glucose is 3 mg/ml.

**catalase is prepared fresh, stock concentration 20 mg/ml in AB buffer.

***glucose oxidase is prepared fresh, stock concentration 25 mg/ml in AB buffer.

8. ATP	Volumes:	1 ml	3 ml	Actual prep:
ATP: always calculate from stock cc.*, e.g.,		3.3 μ l	10 μ l	
GOC		0.9967 ml	2.99 ml	

*ATP stock is 300 mM aqueous solution of ATP, pH 6.8, in 50 μ l aliquots, quick-frozen in liquid nitrogen, and stored at -20°C or, preferably, at -80 °C. Final concentration is adjusted as desired. The numbers are given for 1 mM final ATP concentration.

9. Actin	Volume:	1 ml	5 ml	Actual prep:
TRITC-Ph-F-actin stock		2 μ l	10 μ l ml	
AB/BSA 0.998 ml		4.99 ml		

*Note: dilute actin so that a approximate final dilution of 50,000 of the stored F-actin is obtained. Adjust to a concentration that produces an optimal distribution under the fluorescence microscope.

10. Myosin	Volume:	1 ml	5 ml	Actual prep:
Myosin stock		3 μ l	15 μ l	
0.6 M KCl-AB		0.997 ml	4.985 ml	

*Note: the final concentration of myosin should be about 40 μ g/ml

11. HMM	Volume:	1 ml	5 ml	Actual prep:
HMM stock		0.1 ml	0.5 ml	
AB		0.9 ml	4.5 ml	

*Note: the final concentration of HMM should be about 30-70 μ g/ml

12. S-1	Volume:	1 ml	5 ml	Actual prep:
S-1 stock		0.05 ml	0.25 ml	
AB		0.95 ml	4.75 ml	

*Note: the final concentration of S-1 should be about 50 μ g/ml

13. Additional components of the *in vitro* assay:

a.

b.

c.

Steps:

Check:

1. Prepare flow cells (spin-coat or dip-coat 22x40 mm coverslips with nitrocellulose, and allow to dry overnight. Place two Silastic/Parafilm/glass spacers on the NC surface, and a 18 mm² coverslip on top) _____
2. Pipet aliquots of solutions (150 μ l) in Eppendorf centrifuge tubes _____
3. Degas solutions in desiccator _____
4. Place the flow cell in a 30-ml beaker, and start adding the solutions: _____
(Common sequence:
 - a. 50 μ l myosin (or HMM or S-1) applied from both ends of the flow cell, wait 1 min.
 - b. 100 μ l AB/BSA
 - c. 100 μ l actin, wait 1 min.
 - d. 100 μ l AB/BSA
 - e. 100 μ l GOC)
6. Oil the objective, and place the flow cell on the microscope stage _____
7. Turn on the camera and bring actin filaments in focus. Start recording. _____
8. Add 100 μ l ATP with a micropipet while drawing fluid with filter paper at the other end. _____

Temperature:

Video Recording:

Tape Number:

[illegible]**Notes/Remarks/Modifications:**